


INTRODUCING


ProtekDuo+

With new product enhancements


+ an FDA 510(k) clearance for ECMO




The newest ECMO cannula on the market is one you already know. We've added user-friendly design enhancements to the proven cannula with a 10-year clinical legacy. Experience the next generation with ProtekDuo+.



The **color hemostasis caps** are designed to help visually differentiate the arterial and venous lumens during insertion



The **white printing** on the clamping zone is intended to make it easier to read while there is blood in the cannula



The cannula has **fewer depth insertion markers**, intended to reduce the patient's internal exposure to the ink on the cannula body



THE NEXT GENERATION OF THE PROTEKDUO CANNULA



Ordering Information

To ensure a smooth transition, update your order system with the ProtekDuo+ product numbers.

Size	NEW Product Numbers	Old Product Numbers
31 Fr	5150-5131	5140-5131
29 Fr	5150-4629	5140-4629

Indications for Use:

ECMO: The ProtekDuo+ Veno-Venous Cannula is a single use dual lumen cannula, which provides both venous drainage and reinfusion of blood via the jugular vein, that is indicated for use in adult patients with acute respiratory failure requiring Veno-Venous Extracorporeal Membrane Oxygenation, where other available treatment options have failed and continued clinical deterioration is expected or the risk of death is imminent.

CPB: The ProtekDuo+ Veno-Venous Cannula Set is intended for use as a single cannula for both venous drainage and reinfusion of blood via an internal jugular vein during extracorporeal life support procedures.

Sorin Group Italia S.r.l.
Via Statale 12 Nord, 86
41037 Mirandola (MO) Italy
Tel: +39 0535 29811

Distributed in the USA by:
LivaNova USA
14401 W 65th Way
Arvada, CO 80004 USA



LivaNova
Cardiopulmonary

©2025 LivaNova, Inc. All Rights Reserved. CP-2500048 03/2025